



PHOENIX

Customizable Premium Digital Mobile_x

Turnkey cost-effective OEM Solution.

PHOENIX FULL OEM SOLUTION

- Simple Software/Digital Detector Integration
- High Quality Manufacturing
- Cutting-Edge Technology
- Tight Budget
- Customization
- Regulatory

Bringing your OEM solution with the greatest patient coverage. Phoenix takes advantage of braking energy to recharge batteries while maintaining high image quality with your software/detector integration.

NEW FEATURES TO TRANSFORM DELIVERY CARE



DUAL SIDE COLLIMATOR SHUTTER CONTROL:
for rear/front operator ease.

1 OR 2 TOUCH SCREENS: main 21" monitor with 8,4" secondary tube-head touch-screen (optional).



ENERGY RECOVERY BATTERY TECHNOLOGY:
greatest autonomy to operate without limits.

BENEFITS

- Maximized maneuverability and visibility.
- Up to 25km (15.5mi) and 800 exposures of autonomy.
- Streamline your imaging process.

KEY CUSTOMIZATION

- Free Tube-Head movement with electromagnetic brakes.
- Charging slot for detectors.
- Detector security lock.
- Large storage tray for detectors.
- Status Light and anti-collision sensors.
- Image Preview.
- Up to 50kW. The most powerful generator in the market.
- Remote control for safe exposures.
- Accessories: Barcode reader, DVD, USB, IR-Sensor, RFID Login, tether or wireless connectivity with Hospital network.



EXCELLENCE IN DESIGN IS INCLUDED

- Full View Telescopic Column.
- The smallest footprint.
- Efficient dose management.
- Capacitive handlebar: effortless driving.



CREATING YOUR X-RAY VISION

Your Strategic OEM Partner 



There has never been an easier way to make a New Product Introduction (NPI) and globally distribute your Product. Find your state-of-the-art solution from the comprehensive Sedecal portfolio.

- Complete product process, from design to transfer-to-production, verification & validation, registration & certification and logistic tracking.
- Highest international quality standards, including Quality management system, CE & FDA mark, EMC & electrical safety, X-Ray protection, biocompatibility and sustainability.
- Worldwide registration support, dealing with International organisms and local authorities.
- Heavy investment in R&D, with a cutting-edge engineering and regulatory team, trendy industrial designers and 3D-printing prototyping capability.

Access to all our services and X-ray solutions and SIMPLIFY your NPI on www.sedecal.com



SEDECAL S.A.

C/ Pelaya, 9-13
28110 Algete,
Madrid, SPAIN

+34 916280544
info@sedecal.com
www.sedecal.com

SEDECAL USA, INC

7555 N. Caldwell Ave.
Niles, IL 60714, USA

+001 8473946960
salesoem@sedecalusa.com

BEIJING SEDECAL MEDICAL

Equipment Co., Ltd.
104 Unit, No.1 Building, No.7
Rongchang East Street,
Beijing Economic-technological
Development Area, Beijing, CHINA

+86 10 67806788
info@sedecal.cn
service@sedecal.cn



SEDECAL

www.sedecal.com

